

## Introduction

Human Na<sup>+</sup>,K<sup>+</sup>-ATPase is the target for cardioglycosides such as digitoxin and digoxin which are used in treatment of congestive heart failure and related conditions; thus, it is emerging as an important drug target<sup>1</sup>. The Na<sup>+</sup>, K<sup>+</sup> pump generates electrochemical gradients that are used to drive the coupled transport of many ions and nutrients across the plasma membrane as it actively exports three Na<sup>+</sup> ions with the concomitant import of two K<sup>+</sup> ions, hydrolyzing one ATP molecule in the process. In non-cell-based assays, the activity of Na<sup>+</sup>,K<sup>+</sup>-ATPase has been determined by using purified enzyme preparations to hydrolyze ATP<sup>2</sup>. In cell-based assays, techniques such as patch clamping, fluorescence, H<sup>3</sup>-Ouabain binding, and radio-tracer [Rubidium86], and cold Rubidium flux assay have been used in either recombinant cell lines or in cells other than primary cardiomyocytes<sup>3-5</sup>. However, primary cardiomyocytes cannot be used as they lack the homogeneity, sensitivity and surface binding properties critical to developing cell based assays in a HTS format<sup>6</sup>. In view of the availability of standardized pure cardiac myocytes with functional expression of all essential cardiac ion channels, Aurora Biomed validated Na<sup>+</sup>,K<sup>+</sup>-ATPase in cultured Cor.At® cardiomyocytes derived from transgenic mouse embryonic stem cells.

## Materials & Methods

Rb<sup>+</sup> uptake by cultured Cor.At® cardiomyocytes was carried as follows:

### 1. Culture and maintenance of Cor.At®

**cardiomyocytes:** The cells endogenously expressing Na<sup>+</sup>,K<sup>+</sup>-ATPase were provided frozen in 96 well MTPs, thawed and cultured as per instructions provided by supplier of the cells. On the 4th day of culture, Rb<sup>+</sup> uptake experiment was carried out.

**2. Rb<sup>+</sup> Uptake:** Cells were washed once with 200 uL of Rb<sup>+</sup> Uptake Buffer and Rb<sup>+</sup> uptake was initiated by applying 200 uL of Rb<sup>+</sup> Uptake Buffer.

**a. The Rb<sup>+</sup> uptake profile:** Rb<sup>+</sup> uptake activity of the cells was studied by incubating the cells at specific time intervals in presence of Aurora Biomed's Rb<sup>+</sup> Uptake Buffer.

**b. Dose response:** The cells were incubated in the Rb<sup>+</sup> Uptake Buffer containing appropriate dose of test compound. It was followed by incubation for 15 minutes at room temperature (~22 °C).

**3. Wash:** Residual Rb<sup>+</sup> and compound were removed by four successive washes with 200 uL of SPA-Wash Buffer.


**4. Cell Lysis:** Intracellular samples were obtained by whole cell lysis with the application of 200 uL Lysis Solution.

**5. Analysis:** The level of Rb<sup>+</sup> in the intracellular samples was measured by Aurora Biomed's ICR8000™ using flame atomic absorption spectroscopy.

## Results

**1. Ready to use cells:** The Cor.At® cardiomyocytes were observed to be ready for use for the assay after 3 to 4 days of maintenance.


**2. Expression of Na<sup>+</sup>, K<sup>+</sup>-ATPase isoform:** The Rb<sup>+</sup> uptake profile observed in the assay indicated expression of the this isoform in these cardiomyocytes (Figure 1) showing maximum uptake of Rb<sup>+</sup> in about 80 minutes, with an exponential uptake seen from 0 to 50 minutes. Thus these cells can be employed for 15 minutes of incubation in a HTS-format compound screening assay.



**Figure 1.** Rb<sup>+</sup> uptake profile of Cor.At® cardiomyocytes ATPase isoforms.

### 3. Potency of standard blocker of Na<sup>+</sup>, K<sup>+</sup>-ATPase:

A complete blockage of Na<sup>+</sup>,K<sup>+</sup>-ATPase isoform was observed with 100  $\mu$ M of digitoxin. Digitoxin was determined to have an IC<sub>50</sub> value of 6.82  $\mu$ M (Figure 2). The potency of digitoxin also displays a similar relation between mouse and Chinese hamster Na<sup>+</sup>,K<sup>+</sup>-ATPase isoforms<sup>5</sup>.



**Figure 2.** Cor.At® cardiomyocytes endogenously expressing Na<sup>+</sup>, K<sup>+</sup> -ATPase isoform displaying IC<sub>50</sub> value for digitoxin.

## Conclusion

The data suggest that the Cor.At® cardiomyocytes can be used as a predictive model for identification of modulators of Na<sup>+</sup>, K<sup>+</sup>-ATPase in its original cardiac environment with Rb uptake assay and ICR8000™.

## References

1. De Munari, S. et al. (2003). Structure-based design and synthesis of novel potent Na<sup>+</sup>,K<sup>+</sup>-ATPase inhibitors derived from a 5  $\alpha$ , 14  $\alpha$ androsterane scaffold as positive inotropic compounds. *J Med Chem* 46 (17): 3644-3654.
2. Zahler, R. et al. (1997). Sodium kinetics of Na,K-ATPase  $\alpha$  isoforms in intact transfected HeLa cells. *J Gen Physiol* 110 (2): 201-213.
3. Chen, W. and W. Wu (2002). The asymmetric, rectifier-like I-V curve of the Na/K pump transient currents in frog skeletal muscle fibers. *Bioelectrochemistry* 56 (1-2): 199-202.
4. Geibel, S. et al. (2003). Conformational dynamics of the Na<sup>+</sup>/K<sup>+</sup>-ATPase probed by voltage clamp fluorometry. *Proc Natl Acad Sci U S A* 100 (3):964-969.
5. Gill S, Gill R, Wicks D, Despotovski S & Liang D (2004). Development of an HTS assay for Na<sup>+</sup>, K<sup>+</sup>-ATPase using non-radioactive rubidium ion uptake. *Assay Drug Dev Technol* 2(5):535-42.
6. Automated patch clamp of Cor.At® Cardiomyocytes in the PatchXpress 7000A (MDS-Analytical Technologies) System. Brochure. AxioGeneis AG.