
# STANDARD Q COVID-19 Ag

STANDARD™ Q COVID-19 Ag Test  
PLEASE READ CAREFULLY BEFORE YOU PERFORM THE TEST

REF Q-NCOV-01G

## [ Kit Contents ]

Test device (individually in a foil pouch with desiccant) Extraction buffer tube Filter cap Sterile swab A Sterile swab B (option) Instructions for use

## EXPLANATION AND SUMMARY

### [Introduction]

Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or "2019-nCoV", was discovered because of Wuhan Viral Pneumonia cases in 2019, and was named by the World Health Organization on January 12, 2020, confirming that it can cause colds and the Middle East Respiratory Syndrome (MERS) and more serious diseases such as acute respiratory syndrome (SARS). This kit is helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

### [Intended use]

STANDARD COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of specific antigens to COVID-19 present in human nasopharynx and oropharynx. This test is for *in vitro* professional diagnostic use and intended as an aid to early diagnosis of COVID-19 infection in patient with clinical symptoms with COVID-19 infection. It provides only an initial screening test result. More specific alternative diagnosis methods should be performed in order to obtain the confirmation of COVID-19 infection.

### [Test principle]

STANDARD Q COVID-19 Ag Test has two pre-coated lines, "C" Control line, "T" Test line for the COVID-19 antigen on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimens. Mouse monoclonal anti-COVID-19 IgG antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-COVID-19 IgG antibody conjugated with color particles are used as detectors for COVID-19 antigen device. During the test, COVID-19 antigen in the specimen interact with monoclonal anti-COVID-19 IgG antibody conjugated with color particles making antigen-antibody color particle complex. This complex migrates on the membrane via capillary action until the test line, where it will be captured by the mouse monoclonal anti-COVID-19 IgG antibody recombinant protein. A violet test line would be visible in the result window if COVID-19 antigens are present in the specimen. The intensity of violet test line will vary depending upon the amount COVID-19 antigen present in the specimen. If COVID-19 antigens are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

### [Kit contents]

- Test device (individually in a foil pouch with desiccant)
- Extraction buffer tube
- Filter cap
- Sterile swab A
- Sterile swab B (option)
- Instructions for use

## KIT STORAGE AND STABILITY

Store the kit at room temperature, 2-40°C / 36-104°F, out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

## WARNINGS AND PRECAUTIONS

- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use the extraction buffer tube of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.

## SPECIMEN COLLECTION AND PREPARATION

### [Nasopharyngeal swab]

- To collect a nasopharyngeal swab specimen, insert the sterile swab A into the nostril of the patient, reaching the surface of the posterior nasopharynx.
- Using gentle rotation, push the swab until resistance is met at the level of the turbinate.
- Rotate the swab a few times against the nasopharyngeal wall.
- Remove the swab from the nostril carefully.
- Specimen should be tested as soon as possible after collection.
- If not use of transport media, specimens may be stored at room temperature for up to 24 hours or at 2-8°C / 36-46°F for up to 48 hours in a clean, dry, closed container prior to testing.

### [Throat swab]

- Open the mouth widely.
- Depress the tongue with a tongue based or spoon.
- Insert the sterile swab completely from the mouth into the throat.
- Swab the posterior pharynx, tonsils and other inflamed areas 3 times with moderate force.
- Remove the swab avoid touching to tongue and teeth with the swab.
- Swabs can be stored at room temperature for up to 24 hours or at 2-8°C / 36-46°F for up to 48 hours after collection in a clean and dry plastic tube.

**CAUTION** Minimal dilution of the specimen is recommended, as dilution may result in decreased test sensitivity. **DO NOT USE** a viral transport medium as a specimen. It may cause inaccurate result.

## TEST PROCEDURE

### [Preparation]

- Carefully read instructions for using the STANDARD Q COVID-19 Ag Test.
- Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
- Check the test device and the desiccant pack within the foil pouch.

### [Extraction of specimen]

#### • Nasopharyngeal or Throat swab

- Insert Nasopharyngeal/Throat swab specimen of patient into an extraction buffer tube. Swirl the swab at least five times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Dispose of the used swab in accordance with your biohazard waste disposal protocol.
- Tightly screw the filter cap onto the tube.

### [Analysis of specimen]

- Apply 4 drops of extracted specimen to the specimen well of the test device.
- Read the test result in 30 minutes.

**CAUTION** Positive results should be considered in conjunction with the clinical history and other data available to the physician.

## LIMITATION OF TEST





- The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.
- The test should be used for the detection of COVID-19 antigen in human nasopharyngeal swab or Throat swab samples.
- Neither the quantitative value nor the rate of COVID-19 antigen concentration can be determined by this qualitative test.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- A negative test result may occur if the level of extracted antigen in a sample is below the sensitivity of the test or if a poor quality specimen is obtained.
- For more accuracy of immune status, additional follow-up testing using other laboratory methods is recommended.
- The test result must always be evaluated with other data available to the physician.

## BIBLIOGRAPHY

- Clinical management of severe acute respiratory infection when novel coronavirus(nCoV) infection is suspected. Interim guidance. WHO.2020
- Diagnostic detection of Wuhan coronavirus 2019 by real-time RT-PCR.2020
- Diagnosis and treatment of pneumonia caused by new coronavirus (trial version 4) National Health Commission. 2020

## [Preparation]

- Carefully read instructions for using the STANDARD Q COVID-19 Ag Test.
- Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
- Check the test device and the desiccant pack in the foil pouch.

<Foil pouch>


<Test device>

<Desiccant>

## [Test Procedure]

### Collecting of nasopharyngeal specimen with sterile swab A


- Insert the sterile swab A into the nostril of the patient, reaching the surface of the posterior nasopharynx.
- Swab over the surface of the posterior nasopharynx.
- Withdraw the sterile swab A from the nasal cavity.



**NOTE** Recommendation to test separately for each nasopharyngeal and throat specimen at the same time. In case of discrepancy test result, additional follow-up testing using other laboratory methods is recommended.

### Collecting of throat specimen with Sterile swab B


- Insert the sterile swab B into the throat.
- Swab the posterior pharynx, tonsils and other inflamed areas 3 times with moderate force.
- Withdraw the sterile swab B avoid touching to tongue and teeth with the swab.

- Insert nasopharyngeal specimen swab into an extraction buffer tube. Swirl the swab at least five times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Tightly screw the filter cap onto the tube.
- Apply 4 drops of extracted specimen to the specimen well of the test device.
- Read the test result in 30 minutes.

## [Interpretation of Test Result]

- A colored band will appear in the top section of the result window to show that the test is working properly. This band is control line (C).
  - A colored band will appear in the lower section of the result window. This band is test line of COVID-19 antigen (T).
  - Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.
- \* Positive results should be considered in conjunction with the clinical history and other data available to the physician.**  
**\* The presence of any line no matter how faint the result is considered positive.**

Any inquiries regarding instructions provided should be addressed to: sales@sdbiosensor.com or you can also contact us through www.sdbiosensor.com

Manufactured by SD BIOSENSOR, Inc.  
Head office : C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA  
Manufacturing site : 74, Osongsangmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

Authorized Representative  
MT Promed Consulting GmbH  
Altenhofstrasse 80 66386 St. Ingbert Germany  
Phone : +49 6894 581020, Fax : +49 6894 581021

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**Fabricado por SD BIOSENSOR, Inc.**