





Format: 96 tip
 Tip Volume Range Selection: 10 µL - 1000 µL
 Material: Clear Polypropylene
 Certified RNase-/DNase-free, nonpyrogenic
 Available Options:
 ▶ With or without aerosol-resistant filters
 ▶ Nonsterile or sterilized



Tip Box Format: Non-SBS hanging tip rack
 Material: Number 7 Polycarbonate

Axygen® Robotic Tips

Type	Max Volume	Compatible Head	Tip Box	Packaging	Standard Cat. No.
Non-Filtered	20 µL	8-Channel		Nonsterile Sterile TT-20-C-HTR-S Quantity: 96 tips/rack, 24 racks/case	TT-20-C-HTR
Filtered	10 µL	8-Channel		Sterile Quantity: 96 tips/rack, 24 racks/case	TTF-20-C-HTR-S
Non-Filtered	50 µL	8-Channel		Nonsterile Sterile TT-50-C-HTR-S Quantity: 96 tips/rack, 24 racks/case	TT-50-C-HTR
Filtered	25 µL	8-Channel		Sterile Quantity: 96 tips/rack, 24 racks/case	TTF-50-C-HTR-S
Non-Filtered	200 µL	8-Channel		Nonsterile Sterile TT-200-C-HTR-S Quantity: 96 tips/rack, 24 racks/case	TT-200-C-HTR
Filtered	175 µL	8-Channel		Sterile Quantity: 96 tips/rack, 24 racks/case	TTF-200-C-HTR-S
Non-Filtered	1000 µL	8-Channel		Nonsterile Sterile TT-1000-C-HTR-S Quantity: 96 tips/rack, 16 racks/case	TT-1000-C-HTR
Filtered	1000 µL	8-Channel		Sterile Quantity: 96 tips/rack, 16 racks/case	TTF-1000-C-HTR-S