



# INDICAID<sup>®</sup> OTC COVID-19 Rapid Antigen At-Home Test

## Why Choose INDICAID OTC?

INDICAID OTC COVID-19 Antigen At-Home Tests are a rapid testing device that detects symptomatic and asymptomatic COVID-19 in just 20 minutes. With anterior nasal swabs, specimen collection is quick and comfortable for self-collection. INDICAID OTC is manufactured by Phase Scientific in China and locally distributed by NHD, LLC.

## Features

- Identify current COVID-19 infection status
- FDA Emergency Use Authorized
- Reliable results in 20 minutes
- Simple nasal specimen collection
- Qualitative detection of SARS-CoV-2 antigens
- All-in-one package, including collection swabs



## Kit Contents

- Sold in 12pk or 2pk configuration
- 2/12 test devices
- 2/12 buffer solution vials
- 2/12 sterile nasal swabs
- Instructions for use
- Quick reference guide

## Performance Measures

- 81.7% Sensitivity
- 99.4% Specificity



CALL : 512-651-1977 | EMAIL: [sales@nhd.llc](mailto:sales@nhd.llc)

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | [www.nhd.llc](http://www.nhd.llc)



# INDICAID<sup>®</sup> COVID-19 Rapid Antigen At-Home Test

## Product Specifications and Dimensions

### 2 Pack Configuration

**Kit:** 2 tests / 0.1 lb. / 8" x 1" x 2.6"

**Case:** 252 tests / 126 kits / 16.3 lbs. / 22.4" x 16.9" x 8.7"

**Pallet:** 6,048 tests / 3,024 kits / 389 lbs. / 40" x 40" x 58"

### 12 Pack Configuration

**Kit:** 12 tests / 0.4 lb. / 8.1" x 3.6" x 2.6"

**Case:** 432 tests / 36 kits / 16.3 lbs. / 22.4" x 16.9" x 8.7"

**Pallet:** 10,368 tests / 864 kits / 40" x 40" x 58"

## Attachments and Documentation (click below)

FDA EUA Letter

Instructions for Use

Instructions for Use (Home Test)

Healthcare Provider Fact Sheet



**CALL : 512-651-1977 | EMAIL: [sales@nhd.llc](mailto:sales@nhd.llc)**

5321 Industrial Oaks Blvd., Suite 111 | Austin, TX 78735 | [www.nhd.llc](http://www.nhd.llc)